


3D glasses set(Overglass type)

KM-1010



- 3D glasses for Circular-Polarization Monitor.
- Lightweight body and outstanding visibility.
- Provide protection for your eyes from splash of liquid.
- Exchangeable lens keep the glasses a sanitary condition all the time.

<Specification of KM-1010>

Dimensions: Approx.W160 x H77mm	Operating Atmospheric Pressure: 700~1060hPa
Weight:Approx.23g	Storage/Transportation Temperature: -20~+60℃
Operating Temperature: 0-40℃	Storage/Transportation Humidity: 20~90%
Operating Humidity: 30~85%	Storage/Transportation Atmospheric Pressure: 700~1060hPa



■ 3D Glasses set(Overglass type)
KM-0101(Frame x1,Lens x2)



■ 3D Glasses
AC70(10pcs/set)



■ 26-inch Medical Grade LCD-Color Monitor
MLW-2627C-DC Type-3D

Design and specifications are subject to change without notice.

Ikegami Tsushinki Co.,Ltd.

- Head Office : 5-6-16 Ikegami, Ohta-ku, Tokyo 146-8567, JapanTEL.+81- (0)3-5700-1111 / FAX.+81- (0)3-5700-1137
- Overseas Sales Division : 5-6-16 Ikegami, Ohta-ku, Tokyo 146-8567, JapanTEL.+81- (0)3-5700-4117 / FAX.+81- (0)3-5748-2200
- Beijing Office : Room No.307 Tower B, Jia Tai International Mansion No.41 East Fourth Ring Road, Chaoyang District, Beijing ChinaTEL.+86- (0)10-85710350 / FAX.+86- (0)10-85711804

■ URL : <http://www.ikegami.co.jp/en/>

Ikegami Electronics (U.S.A.),INC.

- Headquarters : 300 Route 17 South,Mahwah,NJ 07430, U.S.ATEL.+1-201-368-9171 / FAX.+1-201-569-1626

■ URL : <http://www.ikegami.com>

Ikegami Electronics (Europe)GmbH

- Headquarters : Ikegami Strasse 1, D-41460 Neuss, GermanyTEL.+49- (0)2131-1230 / FAX.+49- (0)2131-102820
- U.K. Branch : 1 Windmill Business Village, Brooklands Close, Sunbury on Thames, TW16 7DY, EnglandTEL.+44- (0)1932-769700 / FAX.+44- (0)1932-769710

■ URL : <http://www.ikegami.de>

Ikegami Electronics Asia Pacific Pte.Ltd.

- Headquarters : 1 Tampines Central 5, #03-03 CPF Tampines Building, Singapore 529508TEL.+65-6260-8820 / FAX.+65-6260-8896

■ URL : <http://sg.ikegami.co.jp>